

Product datasheet for **TA388828M**


KLHL23 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Kelch-like protein 23 protein (1-200AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q8NBE8
Background:	Cul3-RING ubiquitin ligase complex, protein ubiquitination



[View online »](#)

Product images:

Western blot

All lanes: KLHL23 antibody at 3 μ g/ml

Lane 1: Mouse lung tissue

Lane 2: Human placenta tissue

Lane 3: HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 64 kDa

Observed band size: 64 kDa