

Product datasheet for **TA388850M**

TMEM19 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Transmembrane protein 19 protein (105-217AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q96HH6




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded human lung tissue using [TA388850] at dilution of 1:100



Immunofluorescent analysis of A549 cells using [TA388850] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human gastric cancer using [TA388850] at dilution of 1:100



Western blot
All lanes: TMEM19 antibody at 1 μ g/ml
Lane 1: MCF-7 whole cell lysate
Lane 2: U251 whole cell lysate
Lane 3: Mouse stomach tissue
Lane 4: A549 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 37, 33 kDa
Observed band size: 37 kDa