

## Product datasheet for **TA388852M**


### LLCFC1 Rabbit Polyclonal Antibody

#### Product data:


|                       |                                                                                  |
|-----------------------|----------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                               |
| Applications:         | IHC, WB                                                                          |
| Recommended Dilution: | Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200                           |
| Reactivity:           | Human                                                                            |
| Host:                 | Rabbit                                                                           |
| Isotype:              | IgG                                                                              |
| Clonality:            | Polyclonal                                                                       |
| Immunogen:            | Recombinant Human LLLL and CFNLAS motif-containing protein 1 protein (29-122AA)  |
| Formulation:          | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Concentration:        | lot specific                                                                     |
| Purification:         | >95%, Protein G purified                                                         |
| Conjugation:          | Unconjugated                                                                     |
| Storage:              | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Stability:            | 1 year from dispatch.                                                            |
| Database Link:        | <a href="#">Q96L11</a>                                                           |




[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human liver cancer using [TA388852] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA388852] at dilution of 1:100



Western blot

All lanes: LLCFC1 antibody at 3.5µg/ml + K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 14, 5, 17 kDa

Observed band size: 17 kDa