

Product datasheet for **TA388873M**

C3orf14 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Uncharacterized protein C3orf14 protein (9-128AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9HBI5




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using [TA388873] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using [TA388873] at dilution of 1:100



Western Blot

Positive WB detected in: Rat heart tissue, Mouse heart tissue

All lanes: C3orf14 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa