

Product datasheet for **TA388889**

Accs2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Mouse Acetyl-coenzyme A synthetase, cytoplasmic protein (1-701AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9QXG4
Background:	Activates acetate so that it can be used for lipid synthesis or for energy generation.



[View online »](#)

Product images:



Western Blot

Positive WB detected in: Rat heart tissue, Mouse heart tissue

All lanes: Acss2 antibody at 4 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 79 kDa

Observed band size: 79, 50 kDa