

Product datasheet for **TA388918**

THPO Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Sandwich ELISA: To detect hTPO by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.25 – 1.0 µg/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with ProSci's Polyclonal Anti-Human TPO as a capture antibody, allows the detection of at least 0.2 – 0.4 ng/well of recombinant hTPO. Western Blot To detect hTPO by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 ug/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hTPO is 1.5 - 3.0 ng/lane, under either reducing or non-reducing condition
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hTPO. Anti-Human TPO specific antibody was purified by affinity chromatography and then biotinylated.
Concentration:	lot specific
Purification:	TPO specific antibody was purified by affinity chromatography and then biotinylated
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	P40225



[View online »](#)

Background:

TPO is a lineage-specific growth factor produced in the liver, kidney and skeletal muscle. It stimulates the proliferation and maturation of megakaryocytes, and promotes increased circulating levels of platelets in vivo. TPO signals through the c-mpl receptor, and acts as an important regulator of circulating platelets. Human and murine TPO exhibit cross-species reactivity. The human TPO gene encodes for a 353 amino acid glycoprotein, which contains a 21 amino acid signal sequence, a 15 amino acid erythropoietin-like domain, and a highly glycosylated 179 amino acid C-terminal domain. Recombinant Human TPO is a fully biologically active 174 amino acid polypeptide (18.6 kDa), which contains the erythropoietin-like domain of the full length TPO protein.