

Product datasheet for **TA388964**

DEFB105A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Sandwich ELISA: To detect hBD-5 by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.25 – 1.0 µg/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with ProSci's Polyclonal Anti-Human BD-5 as a capture antibody, allows the detection of at least 0.2 – 0.4 ng/well of recombinant hBD-5. Western Blot To detect hBD-5 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 ug/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hBD-5 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hBD-5. Anti-Human BD-5 specific antibody was purified by affinity chromatography and then biotinylated.
Concentration:	lot specific
Purification:	BD-5 specific antibody was purified by affinity chromatography and then biotinylated
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)

Background:

Defensins (alpha and beta) are cationic peptides with a broad spectrum of antimicrobial activity that comprise an important arm of the innate immune system. The alpha-defensins are distinguished from the beta-defensins by the pairing of their three disulfide bonds. To date, six human beta-defensins have been identified; BD-1, BD-2, BD-3, BD-4, BD-5 and BD-6. beta-defensins are expressed on some leukocytes and at epithelial surfaces. In addition to their direct antimicrobial activities, they can act as chemoattractants towards immature dendritic cells and memory T cells. The beta-defensin proteins are expressed as the C-terminal portion of precursors, and are released by proteolytic cleavage of a signal sequence and, in some cases, a propeptide sequence. beta-defensins contain a six-cysteine motif that forms three intra-molecular disulfide bonds. Recombinant Human BD-5 is a 5.8 kDa protein containing 51 amino acid residues