

## Product datasheet for **TA389074**

### CAV1 Mouse Antibody [Clone ID: M281]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	M281
Applications:	ICC, WB
Recommended Dilution:	<b>WB:</b> 1:500 <b>ICC:</b> 1:50
Reactivity:	Human, Rat, Mouse, Chicken
Host:	Mouse
Isotype:	IgG1
Immunogen:	Clone (M281) was generated from a full length recombinant human Caveolin-1 protein.
Specificity:	This antibody detects a 21 kDa* protein corresponding to the apparent molecular mass of caveolin-1 on SDS-PAGE immunoblots of human A431 cells and rabbit spleen fibroblasts.
Formulation:	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Concentration:	lot specific
Purification:	Protein A Purified
Conjugation:	Unconjugated
Storage:	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Stability:	After date of receipt, stable for at least 1 year at -20°C.
Predicted Protein Size:	21
Database Link:	<a href="#">Q03135</a>



[View online »](#)

**Background:**

Caveolins are the primary structural components of the plasma membrane microdomains, caveolae. Three members of the caveolin family (caveolin-1, -2, and -3) have been identified, and each has distinct expression patterns. Caveolins are involved in diverse biological functions, including vesicular trafficking, cholesterol homeostasis, cell adhesion and apoptosis. Caveolins can interact with various signaling molecules, including G-proteins, receptor tyrosine kinases, PKCs, and Src family kinases. Phosphorylation at Tyr-14 is essential for caveolin association with SH2 or PTB domain-containing adaptor proteins, while phosphorylation at Ser-80 regulates caveolin binding to the ER membrane and entry into the secretory pathway.

**Note:**

Protein G purified tissue culture supernatant.