

## Product datasheet for **TA389103**

### CDK5 Mouse Antibody [Clone ID: DC17+DC34]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	DC17+DC34
Applications:	ICC, IP, WB
Recommended Dilution:	<b>WB:</b> 1:1000 <b>ICC:</b> 1:100
Reactivity:	Human, Rat, Mouse
Host:	Mouse
Isotype:	IgG1
Immunogen:	Clones DC17+DC34 were generated from a recombinant full-length human Cdk5 sequence.
Specificity:	The antibody detects a 31 kDa* band corresponding to Cdk5 on SDS-PAGE immunoblots of human SYF and HeLa cells, and in mouse brain tissue. The antibody does not detect other Cdk family members, such as Cdk1 and Cdk2.
Formulation:	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Concentration:	lot specific
Purification:	Protein A Purified
Conjugation:	Unconjugated
Storage:	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Stability:	After date of receipt, stable for at least 1 year at -20°C.
Predicted Protein Size:	31
Database Link:	<a href="#">Q00535</a>



[View online »](#)

**Background:**

Cyclin-dependent kinases (Cdks) are a family of serine/threonine kinases that require association with regulatory subunits known as cyclins for activation. In addition, post-translational phosphorylation and dephosphorylation events regulate Cdk activity. Phosphorylation of Thr-160 in the T loop by Cdk-activating kinase (CAK) is an obligatory step in kinase activation. By contrast, phosphorylation of the Thr-14 and Tyr-15 residues by the Wee1 family of dual specificity kinases is inhibitory for the Cdks, and dephosphorylation of these residues by the Cdc25 family of phosphatases coincides with Cdk activation. Alternatively, Cdk5 appears to require different mechanisms for activation. This Cdk is activated through association with specific activators, including p35, p39, and p67. Cdk5 is primarily activated in neuronal cells, and only *c-Abl* kinase, rather than Wee family members, have been shown to phosphorylate Tyr-15 to regulate its activity.

**Note:**

Protein G purified tissue culture supernatant.