

Product datasheet for **TA389441**

Biotinylated Anti-UCHL1 antibody(DM214), Rabbit mAb

Product data:

| | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Clone Name: | DM214 |
| Applications: | ELISA, FC |
| Recommended Dilution: | ELISA 1:5000-10000 Flow Cyt 1:100 |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Formulation: | 1XPBS |
| Concentration: | Lot specific |
| Purification: | Purified from cell culture supernatant by affinity chromatography |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature. |
| Stability: | 12 months from date of despatch |
| Background: | The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease.[provided by RefSeq, Sep 2009] |
| Synonyms: | HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1 |



[View online »](#)