

## Product datasheet for **TA389442**

### PE-conjugated Anti-GPRC5D antibody(DM91), Rabbit mAb

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	DM91
Applications:	FC
Recommended Dilution:	Flow Cyt 1:100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	1XPBS
Concentration:	Lot specific
Purification:	Purified from cell culture supernatant by affinity chromatography
Conjugation:	PE
Storage:	+4
Stability:	12 months from date of despatch
Background:	The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.
Synonyms:	GPRC5D



[View online »](#)