

Product datasheet for **TA389595**

Anti-CEACAM6 antibody(DMC685), IgG1 Chimeric mAb

Product data:

| | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | FC |
| Recommended Dilution: | Flow Cyt 1:100 |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Monoclonal |
| Formulation: | 1XPBS |
| Concentration: | Lot specific |
| Purification: | Purified from cell culture supernatant by affinity chromatography |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature. |
| Stability: | 12 months from date of despatch |
| Background: | This gene encodes a protein that belongs to the carcinoembryonic antigen (CEA) family whose members are glycosyl phosphatidyl inositol (GPI) anchored cell surface glycoproteins. Members of this family play a role in cell adhesion and are widely used as tumor markers in serum immunoassay determinations of carcinoma. This gene affects the sensitivity of tumor cells to adenovirus infection. The protein encoded by this gene acts as a receptor for adherent-invasive E. coli adhesion to the surface of ileal epithelial cells in patients with Crohn's disease. This gene is clustered with genes and pseudogenes of the cell adhesion molecules subgroup of the CEA family on chromosome 19. [provided by RefSeq, Apr 2014] |
| Synonyms: | CD66c; CEAL; NCA |



[View online »](#)