

Product datasheet for **TA389599AM**

Anti-GPRC5D (Talquetamab biosimilar) mAb, Biotinylated

Product data:

Product Type:	Primary Antibodies
Applications:	FC
Recommended Dilution:	Flow Cyt 1:10
Reactivity:	Human
Host:	Humanized
Clonality:	Monoclonal
Formulation:	50% Glycerol; PBS, pH 7.4; 0.1% BSA Preservative: 0.1% Procline 300
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Concentration:	Lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Database Link:	Q9NZD1
Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	GPRC5D
Note:	Common Name: JNJ-64407564




Product images:

Figure 1. Flow cytometry analysis with Anti-GPRC5D mAb 25 µg/ml (Talquetamab without CD3) on HEK293 cells transfected with human GPRC5D (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).