

Product datasheet for **TA389602**

Anti-CCR2(plozalizumab biosimilar) mAb

Product data:

| | |
|----------------|--|
| Product Type: | Primary Antibodies |
| Applications: | FC |
| Reactivity: | Human |
| Host: | Humanized |
| Clonality: | Polyclonal |
| Formulation: | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. |
| Concentration: | Lot specific |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature. |
| Stability: | 12 months from date of despatch |
| Database Link: | P41597 |
| Background: | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Synonyms: | CC-CKR-2, CCR2A, CCR2B, CD192, CKR2, CKR2A, CKR2B, CMKBR2, MCP-1-R |



[View online »](#)