

Product datasheet for **TA389604**

Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, FC
Reactivity:	Human
Host:	Humanized
Clonality:	Polyclonal
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Concentration:	Lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.
Stability:	12 months from date of despatch
Database Link:	Q96RD9
Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	BXMAS1; CD307; CD307e; FCRH5; IRTA2; PRO820



[View online »](#)