

Product datasheet for **TA389618**

Rnmt Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB,1:500 - 1:2000 IHC,1:100 - 1:200
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	RNA (guanine-7-) methyltransferase
Database Link:	Entrez Gene 67897 Mouse Q9D0L8
Synonyms:	DKFZp686H1252; hcm1p; hCMT1; hCMT1c; hMet; KIAA0398; MET; OTTHUMP00000162509; RG7MT1



[View online »](#)