

## Product datasheet for **TA389636**

### **FUT11 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fucosyltransferase 11
Database Link:	<a href="#">Entrez Gene 170384 Human Q495W5</a>
Background:	Probable fucosyltransferase.
Synonyms:	Fuc-TXI; FucT-XI; MGC33202; MGC119338; MGC119339



[View online »](#)