

## Product datasheet for **TA389642**

### **YUCCA1 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Arabidopsis
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	<a href="#">Q9SZY8</a>
Background:	Mutant has elevated levels of free IAA in dominant mutant allele; Flavin Monooxygenase-Like Enzyme; Auxin Biosynthesis



[View online »](#)