

Product datasheet for **TA389657**

BSU1 Rabbit pAb Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Arabidopsis
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	Q9LR78
Background:	encodes a serine#Aeithreonine protein phosphatase with an N-terminal Kelch-repeat domain, which is nuclear localized and expressed preferentially in elongating cells. Genetic evidence suggest that this gene plays a redundant role (along with other members of the same gene family) in modulating growth in response to brassinosteroid.



[View online »](#)