

## Product datasheet for **TA389693**

### Aim2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	absent in melanoma 2
Database Link:	<a href="#">Q91VJ1</a>
Background:	Enables double-stranded DNA binding activity and identical protein binding activity. Involved in positive regulation of macromolecule metabolic process; pyroptosis; and regulation of defense response. Acts upstream of or within cellular response to interferon-beta. Located in cytoplasm. Part of AIM2 inflammasome complex. Orthologous to human AIM2 (absent in melanoma 2).
Synonyms:	PYHIN4



[View online »](#)