

Product datasheet for **TA389728**

Ncr1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	natural cytotoxicity triggering receptor 1
Database Link:	Entrez Gene 17086 Mouse Q8C567
Background:	Acts upstream of or within defense response to virus and detection of virus. Located in cell surface. Orthologous to human NCR1 (natural cytotoxicity triggering receptor 1).
Synonyms:	CD335; FLJ99094; hNKp46; LY94; NK-p46; NKP46



[View online »](#)