

Product datasheet for **TA389732**

Thy1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	thymus cell antigen 1, theta
Database Link:	Entrez Gene 21838 Mouse P01831
Background:	This gene encodes a glycoprotein that is anchored to the cell surface of thymocytes, neuronal and other cells through a glycosyl-phosphatidylinositol moiety. A soluble form of the encoded protein has also been detected in serum and cerebrospinal fluid. The encoded protein undergoes further processing to generate the mature protein which mediates cell-cell interactions to trigger downstream signaling pathways.
Synonyms:	CD90; CDw90; FLJ33325



[View online »](#)