

Product datasheet for **TA389733**

Aldh3a1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:1000
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	aldehyde dehydrogenase 3 family, member A1
Database Link:	Entrez Gene 25375 Rat P11883
Background:	Enables 3-chloroallyl aldehyde dehydrogenase activity. Involved in several processes, including response to cAMP; response to glucocorticoid; and response to hypoxia. Located in cytosol. Human ortholog(s) of this gene implicated in arteriosclerosis. Orthologous to human ALDH3A1 (aldehyde dehydrogenase 3 family member A1).
Synonyms:	ALDH3; ALDHIII; MGC10406



[View online »](#)