

Product datasheet for TA389767

Desmoglein 3 (DSG3) Mouse Monoclonal Antibody [Clone ID: 6G2C11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	6G2C11
Applications:	ELISA, FC, IHC, WB
Recommended Dilution:	ELISA (1:10000) Flo (1:200 - 1:400) IHC-P (20 µg/ml) WB (1:200 - 1:2000)
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant fragment of human DSG3 (AA: 55-159) expressed in E. Coli.
Specificity:	Human
Formulation:	PBS, 0.05% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Stability:	12 months from date of despatch
Gene Name:	desmoglein 3
Database Link:	Entrez Gene 1830 Human P32926
Background:	Desmoglein 3 antibody TA389767 is an unconjugated mouse monoclonal antibody to human Desmoglein 3 (DSG3). Validated for ELISA, Flow, IHC and WB. Tested on 20 paraffin-embedded human tissues.
Synonyms:	DSG3, Balding, BAL, Cadherin family member 6, CDHF6, Desmoglein 3, Pemphigus vulgaris antigen, Desmoglein-3, PVA


[View online »](#)