

# **Product datasheet for TA392471**

## **CD147 (BSG) Rabbit Polyclonal Antibody**

### **Product data:**

**Product Type:** Primary Antibodies

Applications: WB

Recommended Dilution: WB: 1:2000~1:5000

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

**Immunogen:** Recombinant protein of human CD147.

**Specificity:** CD147 polyclonal antibody detects endogenous levels of CD147 protein. **Formulation:** Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.

Concentration: 1mg/ml

**Conjugation:** Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Stability: 1 year Predicted Protein Size: ~ 42 kDa

Gene Name: basigin (Ok blood group)

Database Link: Entrez Gene 682 Human

P35613



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



#### Background:

Basigin (EMMPRIN, CD147) is a type I integral membrane receptor protein belonging to the immunoglobulin superfamily. Basigin is a glycosylated protein with four known isoforms, of which isoform 2 is the most abundantly expressed. Multiple functions have been ascribed to Basigin; foremost among these is stimulating the secretion of extracellular matrix metalloproteinases by adjacent fibroblasts, a function which has been implicated in promoting tumor progression. Research studies have shown that Basigin is overexpressed by many tumor cells, and its expression level may correlate with tumor malignancy. A recent study identified the BASIGIN gene as a regulatory target of Slug, suggesting a role for Basigin in the process of epithelial-mesenchymal transition. Basigin has also been identified as a marker for a subset of highly suppressive regulatory T cells, and as an obligate receptor for the malarial parasite Plasmodium falciparum on human erythrocytes.


Synonyms:

5F7; Basigin; BSG; CD147; Collagenase stimulatory factor; EMMPRIN; Extracellular matrix metalloproteinase inducer; HAb18G; Hepatoma-associated antigen; Leukocyte activation antigen M6; OK blood group antigen; TCSF; Tumor cell-derived collagenase stimulatory factor; UNQ6505/PRO21383

Note:


For research use only, not for use in diagnostic procedure.

### **Product images:**



Western blot (WB) analysis of CD147 polyclonal antibody at 1:5000 dilution Lane1:A375 cell membrane lysate(10ug) Lane2:SK-OVCAR3 cell membrane lysate(4ug)





Western blot (WB) analysis of CD147 polyclonal antibody at 1:2000 dilution Lane1:EC9706 cell membrane lysate(25ug) Lane2:A549 cell membrane lysate(25ug)