

Product datasheet for **TA392507S**

IRS1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	WB: 1:10000~1:20000 IF 1:50 - 1:100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to Human IRS-1.
Specificity:	IRS-1 (S1101) polyclonal antibody detects endogenous levels of IRS-1 protein.
Formulation:	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.
Concentration:	1mg/ml
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Stability:	1 year
Predicted Protein Size:	~ 190 kDa
Gene Name:	insulin receptor substrate 1
Database Link:	Entrez Gene 3667 Human P35568
Background:	Insulin receptor substrate 1 (IRS-1) is one of the major substrates of the insulin receptor kinase. IRS-1 contains multiple tyrosine phosphorylation motifs that serve as docking sites for SH2-domain containing proteins that mediate the metabolic and growth-promoting functions of insulin. IRS-1 also contains over 30 potential serine/threonine phosphorylation sites. Ser307 of IRS-1 is phosphorylated by JNK and IKK while Ser789 is phosphorylated by SIK-2, a member of the AMPK family. The PKC and mTOR pathways mediate phosphorylation of IRS-1 at Ser612 and Ser636/639, respectively. Phosphorylation of IRS-1 at Ser1101 is mediated by PKC θ and results in an inhibition of insulin signaling in the cell, suggesting a potential mechanism for insulin resistance in some models of obesity.
Synonyms:	Insulin receptor substrate 1; IRS-1; IRS1




[View online »](#)

Note: For research use only, not for use in diagnostic procedure.

Product images:



Western blot (WB) analysis of IRS-1 (S1101) polyclonal antibody at 1:10000 dilution
Lane1:HepG2 whole cell lysate(40ug) Lane2:MCF-7 whole cell lysate(20ug)



Immunofluorescence analysis of MCF-7 cells using IRS-1 antibody at dilution of 1:100.