

Product datasheet for **TA392682M**

CD11b (ITGAM) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant full length Human CD11B (630aa-1104aa).
Specificity:	CD11B (I1) polyclonal antibody detects endogenous levels of CD11B protein.
Formulation:	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Concentration:	1mg/ml
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Stability:	1 year
Predicted Protein Size:	~ 170 kDa
Gene Name:	integrin subunit alpha M
Database Link:	Entrez Gene 3684 Human P11215
Background:	Integrin α M, also designated complement component receptor-3 α , CD11b (p170), macrophage antigen a polypeptide, cell surface glycoprotein Mac-1 a subunit, MAC1A, MO1A and ITGAM) is a cell adhesion molecule that acts as a receptor for cell surface ligands such as intracellular adhesion molecules (ICAMs) or soluble ligands. Integrins are heterodimeric proteins that contain an α chain and β chain. Integrin α M combines with the Integrin β 2 to form a leukocyte-specific integrin referred to as macrophage receptor 1 (Mac-1), or inactivated-C3b (iC3b) receptor 3 (CR3). Integrin α M/ β 2 is important in the adherence of neutrophils and monocytes to stimulated endothelium, and also in the phagocytosis of complement coated particles.


[View online »](#)

Synonyms: CD11 antigen-like family member B; CD11B; CD11b; Cell surface glycoprotein MAC-1 subunit alpha; CR-3 alpha chain; CR-3 α ; CR3A; CR3 α ; Integrin alpha-M; Integrin α M; ITGAM; Leukocyte adhesion receptor MO1; Neutrophil adherence receptor

Note: For research use only, not for use in diagnostic procedure.

Product images:



Western blot (WB) analysis of CD11B (ITGAM) polyclonal antibody at 1:500 dilution Lane1:THP-1 cell lysate Lane2:Raw264.7 cell lysate Lane3:PC12 cell lysate