

Product datasheet for TA392707

Product datasireet for 1A592707

p16INK4A (CDKN2A) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WE

Recommended Dilution: WB: 1:500~1:1000

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide, corresponding to amino acids N-terminus of Human p16 INK4a.

Specificity: p16 INK4a (A20) polyclonal antibody detects endogenous levels of p16-INK4a protein.

Formulation: Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

Concentration: 1mg/ml

Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Stability: 1 year

Predicted Protein Size: ~ 16 kDa

Gene Name: cyclin-dependent kinase inhibitor 2A

Database Link: Entrez Gene 1029 Human

P42771

Background: The progression of cells through the cell cycle is regulated by a family of proteins designated

cyclin-dependent kinases (Cdks). Sequential activation of individual members of this family and their consequent phosphorylation of critical substrates promotes orderly progression through the cell cycle. Multiple proteins are encoded by the tumor suppressor gene CDKN2A (MTS1/ p16INK4a) via translation of alternate reading frames, resulting in the production of the p19 ARF protein in mice and the p14 ARF protein in humans. p14 ARF induces an increase

in MDM2 and p21 levels and leads to cell cycle arrest in both G1 and G2/M. p14 ARF is negatively regulated by p53 and is known to bind directly to MDM2. CDKN2A also encodes

the mitotic protein p16, which binds to and inhibits the Cdk4/cyclin D complex.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms:

CDK4I; CDKN2; CDKN2A; Cyclin-dependent kinase 4 inhibitor A; Cyclin-dependent kinase inhibitor 2A; MTS-1; MTS1; Multiple tumor suppressor 1; p16-INK4; p16-INK4a; p16INK4A

Note:

For research use only, not for use in diagnostic procedure.

Product images:



Western blot (WB) analysis of p16 INK4a (A20) pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:HEK293T whole cell lysate(40ug) Lane3:PANC1 whole cell lysate(40ug)