

Product datasheet for **TA392807M**

SREB3 (GPR173) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to amino acids 256-310 of Human GPR173.
Specificity:	GPR173 (Y289) polyclonal antibody detects endogenous levels of GPR173 protein.
Formulation:	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Concentration:	1mg/ml
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Stability:	1 year
Predicted Protein Size:	~ 41 kDa
Gene Name:	G protein-coupled receptor 173
Database Link:	Entrez Gene 54328 Human Q9NS66
Background:	G protein-coupled receptors (GPRs or GPCRs), are members of the largest protein family and play a role in many different stimulus-response pathways. G-protein coupled receptors mediate extracellular signals into intracellular signals (G-protein activation). They respond to a great variety of signaling molecules, including hormones, neurotransmitters and other proteins and peptides. GPR173 is also known as super conserved receptor expressed in brain 3 (SREB3). It is an orphan receptor that is expressed primarily in brain and ovary tissues.
Synonyms:	Probable G-protein coupled receptor 173; SREB3
Note:	For research use only, not for use in diagnostic procedure.


[View online »](#)