

## Product datasheet for **TA392975**

### IL2 Polyclonal Antibody [Clone ID: 35C3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	35C3
Applications:	FC
Recommended Dilution:	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Reactivity:	Human
Isotype:	IgG2b, kappa
Clonality:	Polyclonal
Immunogen:	recombinant human IL-2
Specificity:	The mouse monoclonal antibody 35C3 recognizes human interleukin 2 (IL-2; secreted or intracellular).
Formulation:	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide Format: APC Label: Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 2
Database Link:	<a href="#">Entrez Gene 3558 Human P60568</a>
Synonyms:	Interleukin 2



[View online >](#)