

Product datasheet for **TA396621**

Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IP, WB
Recommended Dilution:	WB: User Optimized ELISA: User Optimized
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Turkey serum proteins
Specificity:	This product was prepared from polyspecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-rabbit serum and multiple precipitin arcs against Turkey Serum.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 2.0 mL
Concentration:	90 mg/mL - lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Stability:	Expiration date is one (1) year from date of receipt.
Background:	Anti-Turkey Serum Antibody is specific for the components found in Turkey serum. Turkey Serum provides a broad spectrum of macromolecules, carrier proteins for lipid substances and trace elements, attachment and spreading factors, low molecular weight nutrients, and hormones and growth factors that promote cell growth and health. Turkey serum antibody is suitable for use in veterinary research and biotechnology applications and is ideal for investigators in Cancer and Cell Biology.
Synonyms:	rabbit Anti-turkey whole serum, rabbit anti-turkey serum, anti-serum, detect turkey serum proteins
Note:	Anti-Turkey serum antibody is suitable for use in ELISA, immunoprecipitation, and western blot. Specific conditions for reactivity should be optimized by the end user.



[View online »](#)