

## Product datasheet for **TA397141**

### Chicken IgG Antibody Agarose

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG Antibody Agarose
Applications:	IP
Host:	Goat
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	10 mg/cc - lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti-Chicken IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays.



[View online »](#)