

Product datasheet for **TA397161**

Anti-Goat IgG Antibody Agarose Conjugated

Product data:

Product Type:	Secondary Antibodies
Product Name:	Anti-Goat IgG Antibody Agarose Conjugated
Host:	Rabbit
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	10 mg/mL - lot specific
Conjugation:	Agarose
Storage:	Store Anti-Goat IgG [H&L] (Goat) Antibody Agarose Conjugated vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti-Goat IgG [H&L] (Goat) Antibody Agarose Conjugated is applicable to purification of Goat IgG from sample solution.



[View online »](#)