

Product datasheet for **TA397172**

Guinea Pig IgG Antibody Agarose

Product data:

| | |
|----------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Guinea Pig IgG Antibody Agarose |
| Applications: | IP |
| Host: | Goat |
| Formulation: | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Concentration: | lot specific |
| Conjugation: | Unconjugated |
| Storage: | Store vial at 4° C prior to opening. DO NOT FREEZE. |
| Note: | Anti-Guinea Pig IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays. |



[View online »](#)