

Product datasheet for **TA397184**

Golden Syrian Hamster IgG Antibody Agarose

Product data:

Product Type:	Secondary Antibodies
Product Name:	Golden Syrian Hamster IgG Antibody Agarose
Applications:	IP
Host:	Goat
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti-Golden Syrian Hamster IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays.



[View online »](#)