

Product datasheet for **TA397230**

Human IgG F(c) Antibody Agarose

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG F(c) Antibody Agarose
Applications:	IP
Host:	Rabbit
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti-Human IgG F(c) Agarose Conjugated is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.



[View online »](#)