

## Product datasheet for **TA397245**

### Mouse IgG Secondary Antibody Agarose

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG Secondary Antibody Agarose
Host:	Goat
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	4.4 mg/cc - lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti Mouse IgG Secondary Antibody Agarose can be used for immunoaffinity chromatography of polyclonal or monoclonal mouse IgG antibodies. Once bound, non-specific antibodies are removed by copious washing followed by elution of specific antibody at low pH.



[View online »](#)