

## Product datasheet for **TA397284**

### Rat IgG Antibody Agarose

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG Antibody Agarose
Applications:	IP
Host:	Goat
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	4.4mg/mL - lot specific
Conjugation:	Unconjugated
Storage:	Store Anti-Rat IgG [H&L] (Goat) Antibody Agarose Conjugated vial at 4° C prior to opening. DO NOT FREEZE.
Note:	Anti-Rat IgG [H&L] (Goat) Antibody Agarose Conjugated is applicable to purification of Rat IgG from sample solution.



[View online »](#)