

Product datasheet for **TA397290**

Rat IgM (mu) Antibody Fluorescein Conjugated

Product data:

| | |
|------------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Rat IgM (mu) Antibody Fluorescein Conjugated |
| Applications: | FC, IF |
| Recommended Dilution: | IF: 1:1,000 - 1:5,000 FC: 1:500 - 1:2,500 FLISA: 1:10,000 - 1:50,000 |
| Host: | Rabbit |
| Immunogen: | Rat IgM mu heavy chain |
| Formulation: | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Reconstitution Method: | Restore with deionized water (or equivalent) - Reconstitution Volume: 2.0 mL |
| Concentration: | 10.0 mg/mL - lot specific |
| Conjugation: | FITC |
| Storage: | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| Note: | This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms. |



[View online »](#)