

Product datasheet for **TA397398**

F(ab')₂ Sheep IgG Fab'2 Antibody Peroxidase Conjugated

Product data:

Product Type:	Secondary Antibodies
Product Name:	F(ab') ₂ Sheep IgG Fab'2 Antibody Peroxidase Conjugated
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1:1,000 - 1:10,000 IHC: 1:500 - 1:2,500 ELISA: 1:10,000 - 1:50,000
Host:	Rabbit
Immunogen:	Sheep IgG F(ab') ₂ fragment
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 2.0 mL
Concentration:	lot specific
Conjugation:	HRP
Storage:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note:	F(ab') ₂ Anti-Sheep IgG Fab'2 Peroxidase secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.



[View online »](#)