

## Product datasheet for **TA397667**

### Guinea Pig IgG F(ab')<sub>2</sub> Antibody Peroxidase Conjugated

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Guinea Pig IgG F(ab') <sub>2</sub> Antibody Peroxidase Conjugated
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>WB:</b> 1:1,000 - 1:10,000 <b>IHC:</b> 1:500 - 1:2,500 <b>ELISA:</b> 1:2,000 - 1:6,000
Host:	Goat
Immunogen:	Guinea Pig IgG F(ab') <sub>2</sub> fragment
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 1.0 mL
Concentration:	2.0 mg/mL - lot specific
Conjugation:	HRP
Storage:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note:	Anti-Guinea Pig IgG F(ab') <sub>2</sub> Peroxidase Antibody has been tested by ELISA and is suitable for immunoblotting (western or dot blot), ELISA, immunoelectron microscopy and immunohistochemistry as well as other antibody-based enzymatic assays requiring lot-to-lot consistency.



[View online »](#)