

## Product datasheet for **TA397704**

### Horse IgG F(ab')<sub>2</sub> Antibody Alkaline Phosphatase Conjugated

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Horse IgG F(ab') <sub>2</sub> Antibody Alkaline Phosphatase Conjugated
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>WB:</b> 1:500 - 1:2,500 <b>IHC:</b> 1:200 - 1:1,000 <b>ELISA:</b> 1:2,000 - 1:10,000
Host:	Goat
Immunogen:	Horse IgG F(ab') <sub>2</sub> fragment
Formulation:	0.05 M Tris Chloride, 0.15M Sodium Chloride, 0.001 M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol; pH 8.0
Concentration:	1.0 mg/mL - lot specific
Conjugation:	Alkaline Phosphatase
Storage:	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.
Note:	Anti-Horse IgG F(ab') <sub>2</sub> Alkaline Phosphatase Antibody has been tested by ELISA. This product has been assayed against 1.0 µg of Horse IgG in a standard capture ELISA using pNPP p-nitrophenyl phosphate code # NPP-10 as a substrate for 30 minutes at room temperature. A working dilution of 1:1,000 to 1:4,000 of the reconstitution concentration is suggested for this product.



[View online »](#)