

## Product datasheet for **TA397811**

### Mouse IgG2a Antibody Pre-adsorbed

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a Antibody Pre-adsorbed
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>WB:</b> 1:4,000 - 1:20,000 <b>IHC:</b> 1:2,000 - 1:10,000 <b>ELISA:</b> 1:40,000 - 1:200,000
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2a heavy chain
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	1.1 mg/mL - lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, mix with an equal volume of glycerol, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.
Note:	Anti-Mouse IgG2a secondary antibody has been tested by ELISA and dot blot and is suitable for ELISA, Immunohistochemistry western blotting as well as other anti IgG2a antibody based assays. Mouse IgG2a secondary antibody is available in a variety of formats.



[View online »](#)