

Product datasheet for **TA397835**

Mouse IgG2b Antibody Biotin Conjugated Pre-adsorbed

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b Antibody Biotin Conjugated Pre-adsorbed
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000 ELISA: 1:20,000 - 1:100,000
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2b heavy chain
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 1.0 mL
Concentration:	1.0 mg/mL - lot specific
Conjugation:	Biotin
Storage:	Store second antibody at 4° C prior to restoration. For extended storage aliquot secondary antibody and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note:	Mouse IgG2b biotin conjugated secondary antibodies are available in a variety of formats. Secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.



[View online »](#)