

Product datasheet for **TA397906**

Mouse IgG (gamma 1, 2a, 2b and 3 chain) Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (gamma 1, 2a, 2b and 3 chain) Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000 ELISA: 1:20,000 - 1:100,000
Host:	Rabbit
Immunogen:	highly purified mouse IgG gamma 1, gamma 2a, gamma 2b and gamma 3 proteins
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration:	1.0 mg/ml - lot specific
Conjugation:	Unconjugated
Storage:	Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.
Note:	Anti-Mouse IgG subclass pan reactive Secondary Antibody has been tested by ELISA and is designed for ELISA, Immunohistochemistry, western blotting as well as other anti-IgG antibody based assays. Mouse IgG secondary antibody is available in a variety of formats.



[View online »](#)