

Product datasheet for **TA398124**

F(ab')₂ Dog IgG (H&L) Antibody Phycoerythrin conjugated

Product data:

Product Type:	Secondary Antibodies
Product Name:	F(ab') ₂ Dog IgG (H&L) Antibody Phycoerythrin conjugated
Applications:	FC, IF
Recommended Dilution:	IF: 1:100 - 1:250 FC: 1:100 - 1:250
Host:	Rabbit
Immunogen:	Dog IgG whole molecule
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 1.0 mL
Concentration:	1.0 mg/mL - lot specific
Conjugation:	Phycoerythrin
Storage:	Store vial at 4° C prior to restoration. Restore with deionized water (or equivalent). This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Centrifuge product if not completely clear after standing at room temperature. Do not freeze after reconstitution. Store reagent in the dark. Use subdued lighting during handling and incubation of cells prior to analysis.
Note:	Secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.



[View online »](#)