

## Product datasheet for **TA398177**

### **F(ab')<sub>2</sub> Mouse IgG (H&L) Antibody Biotin Conjugated Pre-Adsorbed**

#### **Product data:**

<b>Product Type:</b>	Secondary Antibodies
<b>Product Name:</b>	F(ab') <sub>2</sub> Mouse IgG (H&L) Antibody Biotin Conjugated Pre-Adsorbed
<b>Applications:</b>	ELISA, IHC, WB
<b>Recommended Dilution:</b>	<b>WB:</b> 1:2,000 - 1:10,000 <b>IHC:</b> 1:1,000 - 1:5,000 <b>ELISA:</b> 1:15,000 - 1:60,000
<b>Host:</b>	Goat
<b>Immunogen:</b>	Mouse IgG whole molecule
<b>Formulation:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Reconstitution Method:</b>	Restore with deionized water (or equivalent) - Reconstitution Volume: 500 µL
<b>Concentration:</b>	1.0 mg/mL - lot specific
<b>Conjugation:</b>	Biotin
<b>Storage:</b>	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Note:</b>	This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms. Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency.



[View online »](#)