

Product datasheet for **TA398185**

F(ab')₂ Mouse IgG (H&L) Antibody Texas Red™ Conjugated Pre-Adsorbed

Product data:

Product Type:	Secondary Antibodies
Product Name:	F(ab') ₂ Mouse IgG (H&L) Antibody Texas Red™ Conjugated Pre-Adsorbed
Applications:	IF
Recommended Dilution:	IF: 1:500-1:2,500
Host:	Donkey
Immunogen:	Mouse IgG whole molecule
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 500 µL
Concentration:	1.0 mg/mL - lot specific
Conjugation:	Texas Red
Storage:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note:	F(ab') ₂ Anti-Mouse IgG Texas Red Antibody has been tested by dot blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.



[View online »](#)