

Product datasheet for **TA398243**

F(ab')₂ Sheep IgG (H&L) Antibody Peroxidase Conjugated

Product data:

Product Type:	Secondary Antibodies
Product Name:	F(ab') ₂ Sheep IgG (H&L) Antibody Peroxidase Conjugated
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1:5,000 - 1:100,000 IHC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:100,000
Host:	Rabbit
Immunogen:	Sheep IgG whole molecule
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reconstitution Method:	Restore with deionized water (or equivalent) - Reconstitution Volume: 500 µL
Concentration:	1.0 mg/mL - lot specific
Conjugation:	HRP
Storage:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note:	F(ab') ₂ Anti-Sheep IgG Peroxidase Conjugated Secondary Antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, as well as other antibody detection methods.



[View online »](#)