

Product datasheet for **TA398290**

Fab Goat IgG (H&L) Antibody Fluorescein Conjugated

Product data:

| | |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Secondary Antibodies |
| Product Name: | Fab Goat IgG (H&L) Antibody Fluorescein Conjugated |
| Applications: | FC, IF |
| Recommended Dilution: | IF: 1:1,000 - 1:5,000 FC: 1:500 - 1:2,500 FLISA: 1:10,000 - 1:50,000 |
| Host: | Rabbit |
| Immunogen: | Goat IgG whole molecule |
| Formulation: | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Reconstitution Method: | Restore with deionized water (or equivalent) - Reconstitution Volume: 500 µL |
| Concentration: | 0.5 mg/mL - lot specific |
| Conjugation: | FITC |
| Storage: | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| Note: | Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. |



[View online »](#)